

	MINISTRY OF RESEARCH, TECHNOLOGY AND HIGHER EDUCATION STATE UNIVERSITY OF JAKARTA FACULTY OF MATHEMATICS AND NATURAL SCIENCES Hasyim Asjari Building, Campus A UNJ Rawamangun Jl. Rawamangun Muka, East Jakarta 13220 Tel/Fax 021-4894909	WORK INSTRUCTION (IK)	
		BIAS INDEX TOOL	
		No. Document	IK ID 01/LF/2022
		Edition	01
		Revision	01
		Is effective	February 2022
Page	1 dari 2		

1. Objective

These instructions are needed as a guide in operating the Refractive Index tool

2. Scope

This work instruction covers the operation of the Refractive Index tool

3. Reference

Refractive Index Practical Book

4. Executor

related PLP.

5. Terms and Definitions

The Refractive Index tool is used to measure reflected or refracted light rays, the result of a light beam hitting the plane between two media.

6. Tool Image



Parts of the Refractive Index Tool

- Power : To turn on the laser
- Knob up/down : To adjust the laser position above the media or below the media
- Glass Tub : To accommodate the liquid media to be measured
- S scale : To determine the direction of the laser towards the media

7. Work instruction

- 7.1 Connect it to a power source
- 7.2 Enter the liquid solution or liquid medium to be measured, make sure the volume is sufficient so that the laser can pass through the medium.
- 7.3 Turn on the tool by pressing the ON button, the indicator light will come on (red and yellow). Set the laser height according to the media to be measured.
- 7.4 Also set the degree deviation that will be used, by moving the refractive index lever.

	MINISTRY OF RESEARCH, TECHNOLOGY AND HIGHER EDUCATION STATE UNIVERSITY OF JAKARTA FACULTY OF MATHEMATICS AND NATURAL SCIENCES Hasyim Asjari Building, Campus A UNJ Rawamangun Jl. Rawamangun Muka, East Jakarta 13220 Tel/Fax 021-4894909	WORK INSTRUCTION (IK)	
		BIAS INDEX TOOL	
		No. Document	IK ID 01/LF/2022
		Edition	01
		Revision	01
		Is effective	February 2022
Page	2 dari 2		

8. Refractive Index Tool Maintenance:

1. Clean the tool from the solution after use.
2. Charge the laser if the light decreases.
3. Store the tool in a cool place with sufficient light and temperature.

9. Endorsement

	Name	Position	Signature	Date
Made by	Nurdi Akbar, S.Pd Muhammad Fajrin S, ST Wulandari Fitriani, M.Pd Muhammad Fajri Z, S.Si Asidiq Saputra, S.Si	Educational Laboratory Institutions		
Checked by	Riser Fahdiran, M.Si	Head of the Physics Laboratory		
Endorsed by	Dr. Widyaningrum Indrasari, M.Sc	Physics Coordination Program		