MODERN PHYSICS PRACTICUM MODULE

Millikan's Oil Drops &
Photoelectric Effect

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UNIVERSITAS NEGERI JAKARTA

PREFACE

In the Name of Allah, the Most Gracious, the Most Merciful, because of His abundance of grace, author can complete this module entitled "MODERN PHYSICS PRACTICE MODULE: MILIKAN'S OIL DROPS & PHOTO-ELECTRIC EFFECT". On this occasion, the author would like to thank all parties who have helped in completing the module that cannot be mentioned one by one. The author really hopes that this module can provide convenience for students in conducting modern physics practicum.

The purpose of making this module is to serve as a companion medium for student practicums in understanding and conducting modern physics practicums, especially in practicum on proprietary oil drops and the photoelectric effect.

Due to the limitations of the author's knowledge, the author realizes that this book has many shortcomings and is far from perfect. Therefore, input from readers means a lot to the author, so that the module can become even better.

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LABORATORY RULES (TATA TERTIB LABORATORIUM)

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- 
- 12) If you experience difficulties or require consultation, you can contact the Laboratory Coordinator.
 - 13) After completing the practicum, tools/materials must be returned to their original place in a complete, clean, and ready-to-use condition. Cleanliness of tools/crockery is the responsibility of the student.
 - 14) Before leaving the laboratory room, the practicum table must be clean and dry, the chairs are arranged in their original places, and the electrical contacts are removed.

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MODULE USE INSTRUCTIONS (INTRUKSI PENGGUNAAN MODUL)

A. Instructions for Students (Petunjuk bagi Mahasiswa)

To obtain maximum learning outcomes, in using this module, the steps taken include:

- a. Carefully read and understand each practicum guide on each topic.
- b. For practical activities, pay attention to the following:
 1. Pay attention to the applicable safety instructions.
 2. Understand each step of the practicum procedure well.
 3. Before carrying out the practicum, prepare the tools and equipment needed.
 4. Ensuring practicum tools function properly.
- c. If you have not mastered the expected level of material, repeat the previous learning activity or ask the lecturer/laboratory assistant concerned.

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B. Petunjuk bagi Laboran (Asisten Laboratorium)

In each practicum activity, the laboratory assistant or laboratory assistant has the role of:

1. Guiding students to carry out practical activities
2. Assisting students in understanding new concepts, practices and answering student questions regarding the practicum process

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3. Helping students to determine and access other additional resources needed for learning.
4. Organize a group learning activities if necessary

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A. Cover (Sampul)

The cover page consists of the title and student identities.

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B. Contents of Report (Isi Laporan)

- 1 Objective
Write the objectives of the practicum in accordance with the experiment being conducted.
- 2 Tools and Materials
Write down the tools and materials used and their quantities.
- 3 theoretical basis
Add theoretical basis from other references/sources such as journals or university-level books (avoid blogs or high

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school books). At least 5 additional references.

4 Work Procedure

Create work procedures in the form of flow charts so that work procedures are not sentences. Use active verbs. Flowcharts are made with charts that have flows that describe the steps or experimental procedures that are made simply, neatly, and clearly.

5 Observation Results

Write down all data for each step taken according to the experiment results. Observation data is written by the sequence of work procedures that have been carried out.

6 Data Processing and Data Analysis

- a. Data processing contains the representation of practicum data, the data that has been obtained is calculated and analyzed into information.
- b. Data analysis contains experimental results and data that has been achieved, calculated according to the theoretical basis, and then compared with the results of existing references.
- c. The results of calculations from experimental data if there is a difference or discrepancy with the theory, it is discussed what causes, obstacles and made based on the theory referred to.
- d. The discussion is written in accordance with the rules of scientific writing. Sentences are written following the rules of good sentence

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writing, which consists of subject, predicate, object, and description. Use various reference sources for comparison

7 Conclusion

The conclusion contains the answer in accordance with the purpose of the experiment written in a simple sentence in the form of bullet points.

8 Bibliography

Write down all references used in accordance with the provisions of writing a minimum of 5 references.

Rules for writing the bibliography:

a. Name

The author's name comes first. The author's last name is written first. This is followed by a comma (,) and then the first and middle names of the author. The title of the author's name does not need to be written.

Example: Dedi Miswardi → Miswardi, D.

If the author consists of two or three people, the author's names are all written with the provision that the first person's name is reversed while the second and third person's names remain. The conjunction "and" is used between the two authors' names.

Example: Dedi Miswardi, Dini Novita, Kartika Sari → Miswardi, D., Dini N., and Kartika S.

If there are more than three people, the first author's name is reversed and the abbreviation "et al" (and

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friends) or et al is added.
Example: Miswardi, D et al.

- b. Year of publication
After the name, include the year of publication of the book used as a reference.
- c. Book title
Write the title of the book in full and write the title in italics.
- d. City and publisher
Write the name of the city first, followed by the publisher's name delimited by a colon (:).
Be sure to use a full stop (.) to delimit the order of the name, year of publication, book title, to the city and publisher's name.

Writing example:

Knight, John F. 2001. Family Medical Care. Bandung: Indonesia Publishing House.

C. Attachment (Lampiran)

The report must be accompanied by an interim report that has been approved by the Practicum Assistant and other supporting attachments if required.

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STATISTICAL METHODS IN DATA PROCESSING (METODE STATISTIK DALAM PENGOLAHAN DATA)

A. Absolute and Relative Uncertainty (Ketidakpastian Mutlak dan Relatif)

The statement $x = (x \pm \Delta x)$ states the absolute uncertainty of the quantity x and states the quality of the measuring instrument used. Relative

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uncertainty can be calculated using the following formula:

$$KSR = \frac{\Delta x}{x} \times 100\%$$

The smaller the relative uncertainty (KSR) value, **the more precise the measurement will be.**

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1. Create Graphics Manually (Membuat Grafik secara Manual)

- Graphics should be made on millimeter block paper
- The graph must be accompanied by a description of the abscissa scale and the ordinate scale
- The graph must be accompanied by a description of the axes
- To connect the points on a graph, the lines must match. If you expect a straight line graph, the line drawn must pass through these points wherever possible (not deviate from the point).
- If you want to draw more than 1 graph in 1 millimeter block, then the points on the graph must be distinguished.

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- The data that has been obtained is copied into Ms. Excel.
- In the Tab Menu Bar section, select Insert, then select Charts.
- In the Charts section there are many chart model choices, select Insert Scatter.

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- d. If a white (empty) box appears, right-click on the mouse or laptop, and select Select Data.
 - e. Click Add, then enter the x-axis and y-axis data (according to the command) by way of the first data block to the end.
 - f. The dots on the graph will appear.
 - g. Click the + sign on the right side of the chart, then check Trendline and Axis Titles.
 - h. Change the chart title and annotate the x-axis and y-axis.
 - i. A straight line will appear between the points.
- d) Jika muncul kotak putih (kosong), klik kanan pada mouse atau laptop, dan pilih Select Data.
 - e) Klik Add, lalu masukkan data sumbu x dan sumbu y (sesuai perintah) dengan cara blok data pertama sampai akhir.
 - f) Titik-titik pada grafik akan muncul.
 - g) Klik tanda + di sisi kanan grafik, lalu centang Trendline dan Axis Titles.
 - h) Ubah judul bagan dan beri anotasi pada sumbu x dan sumbu y
 - i) Garis lurus akan muncul di antara titik-titik tersebut.

C. Linear Regression Method (Metode Regresi Linear (Least Square))

The least squares method is a method that is widely used to see the linear trend of observational data.

The linear relationship between x and y data is:

$$y = mx + c$$

To get the straight line equation above, we can calculate from the following formula:

$$m = \frac{\sum y \sum x^2 - \sum x \sum xy}{n \sum x^2 - (\sum x)^2}$$

$$c = \frac{\sum y - m \sum x}{n}$$

Where y is the ordinate point, x is the abscissa point, m is the gradient, and c is a constant.

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D. Quadratic Regression Method (Metode Regresi Kuadratik)

The quadratic method is the value of the independent variable with a linearly increasing or decreasing shape (occurs parabolically).

General equation of quadratic regression:

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Persamaan umum regresi kudratik:

To find the value of constants a, b, and c, use the following formula:

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
Millikan Oil Drop EXPERIMENT

EKSPERIMEN Tetes Minyak Milikan

In 1909-1913, Robert Andrews Millikan managed to measure the elementary electric charge to unprecedented accuracy and thereby confirmed the quantum nature of the charge. The experiment is based on measuring the quantity of charge carried by charged drops of oil, which are able to rise through the air under the influence of an electric field from a plate capacitor and descend when the field is absent.

Preliminary Stage (Tahap Pendahuluan)

Millikan's Oil Drops
Tetes Minyak Milikan

- Practicum Objectives (Tujuan Praktikum)
- Basic Theory (Teori Dasar)
- Preliminary Assignments (Tugas Pendahuluan)

A. Practicum Objectives (Tujuan Praktikum)

1. Determining the radius of oil (Menentukan jari-jari minyak)
2. Measuring the Elementary Charge (Mengukur Muatan Partikel Elementer)

B. Basic Theory (Teori Dasar)

The principle of the elementary charge of electrons is carried out by observing the motion of oil droplets in an electric field that is formed between two viewing plates as shown in Figure (1.1):

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The motion of oil drops is influenced by gravitational force (F_g), Stoke drag (F_s), Archimedes buoyant force (F_a), and coulomb force (F_c). Observation using a telescope will give an inverted image. Movements without an electric field moving downwards will be observed moving up and motions with an electric field moving up will appear to move downwards.

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1. Motion Without an Electric Field (Gerak Tanpa Medan Listrik)

Note that the direction of the force F_s depends on the direction of the oil drop. The resultant force acting on the oil drop in the absence of an electric field is expressed as follows:

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
Figure 1. 2 1.Motion Without an Electric Field (Gerak Tanpa Medan Listrik)

- b. An object that moves in a fluid will get an upward thrust on Archimedes (F_A):

$$F_A = \frac{4}{3} \pi r^3 \rho_f g \quad (1)$$

- c. The friction of the particles with the fluid will provide a frictional force that is opposite to the direction of motion:

$$F_s = 6\pi \eta r v \quad (2)$$

- d. Particles that have mass will experience a downward gravitational force (F_g) (looks up from the telescope):

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In a balanced state (F_{S1} is a function of v_1), then:

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Keterangan:

r = jari-jari tetes minyak,

ρ_z = density of oil drops,
 ρ_u = mass density of air = 1.29 kg/m³
 v_l = speed of oil drops falling (without an electric field),
 η = viscosity of oil drops in air,
 $\eta = 1.81 \times 10^{-5}$ N s/m²

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If there is an electric field then the oil drops which have a charge of electrons will move/attract upward by the positive plate of the capacitor so the resultant forces acting on the oil drops are:

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Annotation:

v_2 = speed of oil drop moving with the electric field,
 q = the drop charge of oil,
 V_{AB} = voltage difference between the capacitor plates and
 d = distance traveled by the oil drop moving up.

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The oil drop charge is a multiple of the base charge value e and is expressed by the following equation:

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$$q = ne \tag{7}$$

With $e = 1,6 \times 10^{-19}$ C, n (integer) = 1, 2, 3, 4

When the oil drops move in the air medium, evaporation can occur due to friction increasing air temperature which affects the viscosity value. The correction factor (Cunningham factor) viscosity value is expressed by the following equation [8]:

$$\eta_c(T) = \eta_0(T) \left[1 + \frac{b}{rP} \right]^{-1} \tag{8}$$

where η_c is the drop-corrected viscosity of oil, η_0 is the viscosity of air as a function of T, b (constant) = 6.17×10^{-4} and P is air pressure.

Correction of the viscosity value results in a correction in determining the value of the oil drop charge which is expressed by the following equation [9]:

$$q_c = q \left[1 + \frac{b}{rP} \right]^{-3/2} \tag{9}$$

3. Observational Method (Metode Pengamatan)

a. *Floating method* (Metode mengambang)

The voltage U which keeps the drop floating still is determined along with the velocity v_2 at which the drop descends after the voltage is turned off

$$v_2 = \frac{x}{t_2} = \frac{S}{V \cdot t_2} \tag{10}$$

Annotation:

- t_2 : Time to descend,
- S: Separation of scale positions,
- V: Objective magnification (2x)

The radius r_0 and charge q_0 of the oil drop are determined from the equilibrium level of electrical force, buoyancy in air, Stokes' friction in air and the force of gravity:

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Annotation:

η : Viscosity of air,
 ρ_2 : Density of oil,
 ρ_1 : Density of air
 g : Acceleration due to gravity
 d : Separation of capacitor plates
 (3 mm)

Keterangan:

η : Viskositas udara,
 ρ_2 : kerapatan minyak
 ρ_1 : kerapatan udara
 g : gaya gravitasi
 d : Pemisahan plat kapasitor (3 mm)

b. *Rising/falling method* (Metode naik/turun)

For a selected value of voltage U , the velocity with which a drop rises v_1 and the velocity with which it descends v_2 are both determined

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Annotation:

η : Viscosity of air,
 ρ_2 : Density of oil,
 ρ_1 : Density of air
 g : Acceleration due to gravity
 d : Separation of capacitor plates (3 mm)

Keterangan:

η : Viskositas udara,
 ρ_2 : kerapatan minyak
 ρ_1 : kerapatan udara
 g : gaya gravitasi
 d : Pemisahan plat kapasitor (3 mm)

4. Viscosity and Stokes' Law / Viskositas dan Hukum Stokes

Density of oil at temperature 15°C
and 25°C:

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1. Adjustable feet
2. Measuring microscope on stand
3. Light intensity adjustment
4. Voltage switch U
5. Timer switch t
6. Experiment chamber
7. Display and control unit
8. Oil atomiser
9. Switch for changing plate capacitor sign
10. Co-axial power socket for plug-in power supply
11. Voltage adjustment
12. Rubber bulb
13. Focussing knobs

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
Figure 1.4 Experiment chamber: 1 Cover, 2 Hole for oil atomiser, 3 Top capacitor plate, 4 Observation window, 5 Vertical adjustment knob for head of microscope, 6 Stand rod for measuring microscope

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Note: The goal is to generate a small number of individual oil droplets, not an extensive collection of droplets as you only need to select a single droplet for observation. The matched droplet will appear as a brightly lit dot on the meter of the focusing microscope. Spraying oil more than once resulted in too many droplets in the chamber, especially in the area

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between the microscope and the focus. A large number of drops will interfere with the drop protection that is in the direction. If too much oil gets into the chamber, it must be cleaned. If no droplets appear on the microscope, even after several sprayings of oil, the openings in the top plate may be clogged. Then it needs to be cleaned too.

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I. *Floating method* (Metode mengambang)

Practicum steps taken:

1. Reset the time to zero, If the time has been saved from before, by pressing the reset screen.
2. Determines the upper and lower threshold distances in the metering lens
3. Choose the voltage polarity U , and watch the movement of the oil drops.

If the position of the voltage polarity is above, it means that the capacitor plates are positive above and negative below. Drops of positive particles will be attracted by a negative field. However, if the motion of the particles remaining downward means that the particles are charged with a direct voltage, then it is necessary to change the polarity of the voltage downwards for the polarity to change.

4. Turns on the voltage and timer switch to ON. and ignore the time at t_1

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II. *Rising/falling method* (Metode naik/turun)

Practicum steps taken:

- 1 Reset the time to zero, If the time has been saved from before, by pressing the reset screen.
- 2 Determines the upper and lower threshold distances in the metering lens
- 3 Choose the voltage polarity U , and watch the movement of the oil drops.

If the position of the voltage polarity is above, it means that the capacitor plates are positive above and negative below. Drops of positive particles will be attracted by a negative field. However, if the motion of the particles remaining downward means that the particles are charged with a direct voltage,

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E. Data Analysis of Practicum Results (Analisis Data Hasil Praktikum)

1. Data Processing / Pengolahan Data

(a) Calculate the relative uncertainty of the voltage in the experiment with the floating and up/down method and write down the measurement results! Record the results in table 1.3.

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(b) Calculate the relative uncertainty of the voltage in the experiment with the floating and up/down method and write down the measurement results! Record the results in table 1.4.

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(c) Calculate the Velocity ($v = \frac{s}{t}$) of Oil Drops, record in table 1.5

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(d) Calculate the radius of the oil drop using eq. 5, and record it in the following table 1.5

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
chart 1. 1 The relationship between voltage and time fall/rising method (Hubungan tegangan terhadap waktu metode naik/turun)

- (b) Based on table 1.2, draw a graph of Time (t) against the voltage (U)!
- (b) Berdasarkan tabel 1.2, buatlah Grafik waktu (t) terhadap tegangan (U)!


chart 1. 2 The relationship between the radius of the oil drop (r) and the velocity (v) (Hubungan Jari-jari tetes minyak (r) terhadap kecepatan(v))

(c) Based on table 1.3, draw a graph of the radius of the oil drop (r) against the velocity (v)!

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Answer/Jawab:

- (c) Based on graphs 1 and 2, what is the relationship between time and voltage? (Berdasarkan grafik 1 dan 2, bagaimanakah hubungan waktu dan tegangan?)

Answer/Jawab:

- (d) Based on graph 3, what is the relationship between the radius of the oil drop (r) and the velocity (v)? (Berdasarkan grafik 3, bagaimanakah hubungan Jari-jari tetes minyak (r) terhadap kecepatan(v)?)

Answer/Jawab:

4. Kesimpulan

- (a) Draw conclusions about the relationship between stress and time, viscosity values and particle velocity based on the analysis of the experimental results! (Buatlah kesimpulan hubungan antara tegangan terhadap waktu, nilai viskositas dan kecepatan partikel berdasarkan analisis hasil percobaan yang dilakukan!)

Answer/Jawab:

- (b) Make a conclusion about the value of the radius of the load and Elementary charge based on the calculation of the data processing! (Buatlah kesimpulan nilai jari-jari muatan dan Muatan Elementer berdasarkan perhitungan pada pengolahan data!)

Answer/Jawab:

REFERENCE

Alvensleben, L.V. Phywe experimental literature physics: Elementary charge and Millikan experiment. LEP 5.1.01

Alvensleben, L.V. Phywe experimental literature physics: Elementary charge dan Milikan experiment. LEP 5.1.01

Gandarilla, Iran. 2022. “Millikan Oil Drop Experiment” *Digital Showcase*. Lynchburg: University of Lynchburg. hlm 3

Leybold instruction sheet 559 41/42. Milikan apparatus power supply.

Melissinos, A.C. (1966). Experiment in modern physics. New York: Academic press.

Milikan, R.A. (1913). On the elementary electrical charge and the Avogadro constant. Physical Review, Vol. II, No. 2, pp. 109-143.

Raymond A. Serway, Clement J. Moses and Curt A Moyer. (2005). Modern physics, 3 rd edition, Belmont: Thomson learning, Inc.

Sharma, Shapna dan Ahluwalia, P K. 2018. “Can virtual labs become a new normal? A case study of Millikan’s oil drop experiment” *European Journal of Physics XXXVIII*

Photoelectric Effect **EXPERIMENT**

EKSPERIMEN *Efek Fotolistrik*

Light of known frequency passes through a ring-shaped anode to collide with a cathode, where it causes electrons to be released due to the photoelectric effect. The energy of the electrons can be determined by applying a decelerating voltage, which compensates for the flow of electrons towards the anode until no electrons are flowing. The energy of the electrons is therefore similarly independent of intensity.

By obtaining the cut-off voltages for light of varying frequency, it is possible to calculate Planck's constant.

Preliminary Stage (Tahap Pendahuluan)

Photoelectric Effect
Efek Fotolistrik

- Practicum Objectives (Tujuan Praktikum)
- Basic Theory (Teori Dasar)
- Preliminary Assignments (Tugas Pendahuluan)

A. Practicum Objectives / Tujuan Praktikum

- | | |
|--|---|
| <ol style="list-style-type: none"> 1. Determine the value of the Planck constant 2. Calculate the wavelength of light 3. Verification of the inverse square law of radiation using a photoelectric cell | <ol style="list-style-type: none"> 1) Menentukan nilai tetapan Planck 2) Mempelajari efek fotolistrik 3) Verifikasi hukum kuadrat terbalik radiasi menggunakan sel fotolistrik |
|--|---|

B. Basic Theory/Teori Dasar

Einstein's (1905) postulate was that light consists of a flux of particles, called photons, whose energy is proportional to frequency through the relationship:

$$E = h\nu$$

The proportionality factor h is known as Planck's constant, and is considered a natural constant. In this concept, each photo electron is liberated by light (photon) and exits the atom with the maximum kinetic energy expressed by:

$$dK = h\nu - W_0$$

work function whose value depends on the type of metal

The next photo electrons move to the anode and will form a photo electron current. If the negative potential of the anode continues to increase, the photoelectron current will decrease until it becomes zero. The voltage at which the photoelectron current is zero is called the limit voltage. The relationship between the kinetic energy is:

$$eU_0 = h\nu - W_0$$

$$eU_0 = h\nu$$

With $e = 1.6 \times 10^{-19} \text{C}$. In this state all photo electrons even with maximum kinetic energy cannot reach the anode so that the quantity is no longer a work function.

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The relationship between the frequency of light and the wavelength of light is expressed by the equation:

$$v = c/\lambda$$

Where c is the speed of light (3×10^8 m/s). When the wavelength of light is changed so that the frequency of light increases by Δv , the energy of the photo electrons will increase by $h\Delta v$ so that the limit voltage can be increased by ΔU_0 . Thus the relationship is obtained:

$$\frac{\Delta U_0}{\Delta v} = \frac{h}{e}$$

Einstein's theory predicts that the frequency of the incident light varies, and the "stopping" potential, V , plotted as a function of frequency gives the slope of the line h/e (see Figure 2.1).


Figure 2. 1 The slope of the line h/e
Kemiringan Garis h/e

Verifying the inverse square law of radiation using a photoelectric cell.

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If L is the light intensity of the electric lamp and E is the illuminance (illumination intensity) at point " r " of it, then according to the inverse square law.

$$E = \frac{L}{r^2}$$

If light is allowed to face the cathode of a photoelectric cell, the photoelectric current (I) will be proportional to E .

$$E = \frac{L}{r^2} = K \cdot I$$

Therefore the graph between $\frac{1}{r^2}$ and I is a straight line, which verifies the inverse square law of radiation.

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C. Video of the Concept Application of the photoelectric effect / Video Penerapan Konsep Efek Fotolistrik

To find out the application of the photoelectric effect, please scan the following QR-Code and observe the video provided. (Untuk mengetahui penerapan efek fotolistrik, silakan pindai QR-Code berikut ini dan amati video yang disediakan).

Figure 2. 2 QR Code Application of the photoelectric effect (Kode QR Aplikasi Efek Fotolistrik)

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
Figure 2. 3 the photoelectric effect tools
Alat efek fotolistrik

Annotation:

1. Light Source
2. The arena of movement of light
3. Ruler scale
4. Drawtube (Where the color filter is placed)
5. Covers
6. Focus Lens
7. Vacuum phototube
8. Base for holding the phototube
9. Digital Meters

Keterangan Gambar:

1. Sumber Cahaya
2. Arena pergerakan cahaya
3. Skala penggaris
4. Drawtube (Tempat filter warna ditempatkan)
5. Cover
6. Lensa Fokus
7. Phototube hampa udara
8. Dasar untuk memegang phototube

10. Display mode
11. Current range setting
12. Setting the light intensity
13. Color filters
14. Voltage regulation
15. Lens cover
16. Voltage direction setting
17. Power
18. Power Indicator

9. Digital Meter
10. Display mode
11. Pengaturan rentang arus
12. Pengaturan intensitas cahaya
13. Filter warna
14. Pengaturan tegangan
15. Lensa penutup
16. Pengaturan arah tegangan
17. Power
18. Power Indicator

3. Practicum Steps/Langkah Kerja

Scan the QR-Code on the side to display the practicum video (Pindai QR-Code di samping untuk menampilkan video praktikum).

a. Tool Calibration/Kalibrasi Alat

1. Make sure the appliance is plugged into a power supply
 2. Adjust the position of the lamp with a distance of 35 cm
 3. Make sure the lens is closed
 4. Turn on the photoelectric effect tool
 5. Adjust the light intensity with medium conditions
 6. Check the current flowing in the closed lens condition, make sure no current flows or rated current $<0.003\mu\text{A}$
 7. Open the lens cover and check for current flowing without a filter
 8. Record the value of the flowing current
- 1) Pastikan kabel alat telah terhubung dengan stop
 - 2) Aturlah posisi lampu dengan jarak 35 cm
 - 3) Pastikan lensa dalam keadaan tertutup
 - 4) Nyalakan alat efek fotolistrik
 - 5) Aturlah intensitas cahaya dengan kondisi medium
 - 6) Ceklah arus yang mengalir pada kondisi lensa tertutup, pastikan tidak ada arus yang mengalir atau nilai arus $<0,003\mu\text{A}$
 - 7) Bukalah penutup lensa dan cek arus yang mengalir tanpa filter
 - 8) Catatlah nilai arus yang mengalir

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2. Grafik Hubungan Antar Variabel

- a. Based on table 2.3, make a graph of frequency (ν) against stopping potential (ν_s)

- a) Berdasarkan tabel 2.3, buatlah Grafik frekuensi (ν) terhadap stopping potential (ν_s)


Chart 2. 1 Relationship of Frequency (ν) with stopping potential (U) (Hubungan Frekuensi (ν) dengan tegangan henti (U))

- b. Based on table 2.5, make an inverse graph of the distance squared ($\frac{1}{r^2}$) to the electric current (I)

- b) Berdasarkan tabel 2.5, buatlah Grafik terbalik jarak kuadrat ($\frac{1}{r^2}$) terhadap arus listrik (I)


Chart 2. 2 The relationship of the squared distance ($\frac{1}{r^2}$) to the electric current (I) (Hubungan jarak kuadrat ($\frac{1}{r^2}$) terhadap arus listrik (I))

3. Graph Analysis (Analisis Grafik)

- (a) Based on graphs 1, is frequency linear with respect to stopping potential? What caused this to happen? (Berdasarkan grafik 1, apakah frekuensi linear terhadap tegangan henti? Apa yang menyebabkan terjadinya hal tersebut?)

Answer/Jawab:

- (b) Based on graph 2, what is the distance (r) to the electric current (I)? What caused this to happen? (Berdasarkan grafik 3, apakah jarak (r) terhadap arus listrik (I)? Apa yang menyebabkan terjadinya hal tersebut?)

Answer/Jawab:

4. Conclusion (Kesimpulan)

- (a) Make a conclusion about the relationship between voltage, current and distance based on the analysis of the experimental results! (Buatlah kesimpulan hubungan antara tegangan kuat arus dan jarak berdasarkan analisis hasil percobaan yang dilakukan!)

Answer/Jawab:

- (b) Make a conclusion about the Planck constant value obtained based on calculations in data processing! (Buatlah kesimpulan nilai konstanta planck yang didapatkan berdasarkan perhitungan pada pengolahan data!)

Answer/Jawab:

REFERENCE

Leybold Physics leaflet. Determine Planck’s constant. Leybold Didactic GMBH.
Raymond A. Serway, Clement J. Moses and Curt A Moyer. (2005). Modern physics, 3 rd edition, Belmont: Thomson learning, Inc.
Suliyana dkk. Panduan Praktikum Fisika Modern. Semarang:JDS